

ICX-Systemoverview 6.0

... reloaded ...



59,-€*
per ICX-implant
all lengths
all diameters
*plus VAT

SILVERSPONSOR
of DGI

BRONZESPONSOR
of EAO

ICX-IMPLANTS	
Overview ICX-Implants	3
ICX-HEALING CAPS	
Overview ICX-Cover Screw- and Healing Caps	4
Overview ICX-Healing Caps Individual & Customizable	5
ICX- IMPRESSION POSTS	
Overview ICX-Titanium Impression Posts	6
ICX-LAB ANALOGS	
Overview ICX-Lab Analogs.....	7
ICX-ABUTMENTS	
Overview ICX-Titanium-Abutments	8
Overview ICX-Massive-Abutments.....	10
Overview ICX-Universal-Abutments	10
Overview ICX-Bar-, Gold- and Plastic-Abutments	11
Overview ICX-plus-Abutments.....	12
Overview ICX-CAD/CAM-Adhesive- and Scanabutments.....	12
Overview ICX-Ceramic-Abutments.....	13
Overview ICX-Maximus-Abutments, twopiece	14
Overview ICX-Maximus-Abutments, onepiece	14
Overview Locator™ for the ICX-templant®	15
ICX-t-bona	
Overview ICX-t-bona Ball-Head-System.....	16
Overview ICX-mini Ø 2.9mm/Prosthetics.....	17
OCC-SYSTEM	
Overview OCC-Product-Line =occlusally screw-retained implant prosthetics.....	18-19
ICX-DRILLS	
Overview ICX-Parallel- and Stop-Drills.....	20
Overview ICX-Drillstops, Titanium	22
ICX-Drillstops-Box.....	23
ICX-ACCESSOIRES	
Overview ICX-templant®-Accessoires	24-25
ICX-Surgery-Boxes	26

Overview: ICX-Implants


ICX-mini Implant white Ø 2.9mm

per piece: € 59,-/pack of 6 only: € 354,-

Article-Nr.

C-001-290080	ICX-mini Ø 2.9 x 8mm (pack of 6 only)
C-001-290100	ICX-mini Ø 2.9 x 10mm (pack of 6 only)
C-001-290125	ICX-mini Ø 2.9 x 12.5mm (pack of 6 only)

Note: Only use C-003-000004 ICX-mini Cover cap!




ICX-plus Implant yellow, Ø 3.45mm

per piece: € 59,-/pack of 6 only: € 354,-

Article-Nr.

C-001-345065	ICX-plus Ø 3.45 x 6.5mm (pack of 6 only)
C-001-345080	ICX-plus Ø 3.45 x 8mm (pack of 6 only)
C-001-345100	ICX-plus Ø 3.45 x 10mm (pack of 6 only)
C-001-345125	ICX-plus Ø 3.45 x 12.5mm (pack of 6 only)




ICX-Implant red, Ø 3.75mm

per piece: € 59,-/pack of 6 only: € 354,-

Article-Nr.

C-001-390080	ICX Ø 3.75 x 8mm (pack of 6 only)
C-001-390100	ICX Ø 3.75 x 10mm (pack of 6 only)
C-001-390125	ICX Ø 3.75 x 12.5mm (pack of 6 only)
C-001-390150	ICX Ø 3.75 x 15mm (pack of 6 only)




ICX-Implant green, Ø 4.1mm

per piece: € 59,-/pack of 6 only: € 354,-

Article-Nr.

C-001-425080	ICX Ø 4.1 x 8mm (pack of 6 only)
C-001-425100	ICX Ø 4.1 x 10mm (pack of 6 only)
C-001-425125	ICX Ø 4.1 x 12.5mm (pack of 6 only)
C-001-425150	ICX Ø 4.1 x 15mm (pack of 6 only)




ICX-Implant blue, Ø 4.8mm

per piece: € 59,-/pack of 6 only: € 354,-

Article-Nr.

C-001-480080	ICX Ø 4.8 x 8mm (pack of 6 only)
C-001-480100	ICX Ø 4.8 x 10mm (pack of 6 only)
C-001-480125	ICX Ø 4.8 x 12.5mm (pack of 6 only)
C-001-480150	ICX Ø 4.8 x 15mm (pack of 6 only)



*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Overview: ICX-Cover Screw- and Healing Caps

Article-No.	Article	Dimensions	Material	Price Euro	
C-003-000001	ICX-Cover Screw	-	Titanium	16,50*	
C-003-000002	ICX-Cover Screw, sterile, pack of 6 pieces only	-	Titanium	99,-*	
C-004-444720	ICX-Healing Cap Titanium, bottle size Ø 3.9/4.7mm	2 mm Gingiva Height 3.5mm H	Titanium	29,-*	
C-004-444730	ICX-Healing Cap Titanium, bottle size Ø 3.9/4.7mm	3 mm Gingiva Height 4.5mm H	Titanium	29,-*	
C-004-444740	ICX-Healing Cap Titanium, bottle size Ø 3.9/4.7mm	4 mm Gingiva Height 5.5mm H	Titanium	29,-*	
C-004-444750	ICX-Healing Cap Titanium, bottle size Ø 3.9/4.7mm	5 mm Gingiva Height 6.5mm H	Titanium	29,-*	
C-004-004510	ICX-Healing Cap Ø 4.5mm conical	1 mm Gingiva Height	Titanium	29,-*	
C-004-004520	ICX-Healing Cap Ø 4.5mm conical	2 mm Gingiva Height	Titanium	29,-*	
C-004-004530	ICX-Healing Cap Ø 4.5mm conical	3 mm Gingiva Height	Titanium	29,-*	
C-004-004540	ICX-Healing Cap Ø 4.5mm conical	4 mm Gingiva Height	Titanium	29,-*	
C-004-004550	ICX-Healing Cap Ø 4.5mm conical	5 mm Gingiva Height	Titanium	29,-*	
C-004-004560	ICX-Healing Cap Ø 4.5mm conical	6 mm Gingiva Height	Titanium	29,-*	
C-004-006020	ICX-Healing Cap Ø 6.0mm conical	2 mm Gingiva Height	Titanium	29,-*	
C-004-006030	ICX-Healing Cap Ø 6.0mm conical	3 mm Gingiva Height	Titanium	29,-*	
C-004-006040	ICX-Healing Cap Ø 6.0mm conical	4 mm Gingiva Height	Titanium	29,-*	
C-004-006050	ICX-Healing Cap Ø 6.0mm conical	5 mm Gingiva Height	Titanium	29,-*	
C-004-006060	ICX-Healing Cap Ø 6.0mm conical	6 mm Gingiva Height	Titanium	29,-*	

Note: Please tighten all ICX-cover screws and healing caps by hand only. A torque of 5 to 10Ncm should not be exceeded. There is a risk that cone forces may lock it permanently, if screw with 20 Ncm or more

Note: H = Height from implant shoulder

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Overview: ICX-Healing Caps Individual & Customizable

Article-No.	Article	Dimensions	Material	Price Euro	
C-004-504720	ICX-Healing Cap, individual with hex, Ø 4.7mm, H 2, Peek	2mm H 6mm K	Peek, sterilizable	29,-*	
C-004-504740	ICX-Healing Cap, individual with hex, Ø 4.7mm, H 4, Peek	4mm H 6mm K	Peek, sterilizable	29,-*	
C-004-514720	ICX-Healing Cap, individual non hex, Ø 4.7mm, H 2, Peek	2mm H 6mm K	Peek, sterilizable	29,-*	
C-004-514740	ICX-Healing Cap, individual non hex, Ø 4.7mm, H 4, Peek	4mm H 6mm K	Peek, sterilizable	29,-*	
C-004-524706	ICX-Healing Cap, customizable with hex, Ø 5mm, H 6, Peek	6mm H	Peek, sterilizable	29,-*	
C-004-524708	ICX-Healing Cap, customizable with hex, Ø 5mm H 8, Peek	8mm H	Peek, sterilizable	29,-*	
C-004-524710	ICX-Healing Cap, customizable with hex, Ø 5mm, H 10, Peek	10mm H	Peek, sterilizable	29,-*	
C-004-526006	ICX-Healing Cap, customizable with hex, Ø 7mm, H 6, Peek	6mm H	Peek, sterilizable	29,-*	
C-004-526008	ICX-Healing Cap, customizable with hex, Ø 7mm, H 8, Peek	8mm H	Peek, sterilizable	29,-*	
C-004-526010	ICX-Healing Cap, customizable with hex, Ø 7mm, H 10, Peek	10mm H	Peek, sterilizable	29,-*	





*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)



All articles on this page: Connection Screw included

Note: H = Height from implant shoulder
K = Length of the stump starting from preparation margin




Overview: ICX-Titanium Impression Posts

Article-No.		Article	Dimensions	Material	Price Euro
C-005-020002		<p>ICX-Titanium Impression Post, closed impression, 3 pieces consisting of: closed Impression Coping Connection Screw Impression Cap</p> <p>Note: this article replaces the older article with same article number</p>	11.5mm H incl. cap	Titanium	39,-*
C-005-040002		ICX-Impression Cap, blue, pack of 10 pieces only		Plastic	59,-*
C-005-030001		<p>ICX-Titanium Impression Post, open impression, 2 pieces contains: open Impression Coping Connection Screw</p>	4mm H	Titanium	39,-*
C-005-030002		<p>ICX-Titanium Impression Post, open impression, 2 pieces, Ø 4.5mm contains: open Impression Coping Connection Screw</p>	11mm H	Titanium	39,-*

Note: H = Height from implant shoulder

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

















Overview: ICX-Lab Analogs

C-006-010001		ICX-Lab Analog	15mm Length	Titanium	19,-*
C-006-010002		ICX-plus-Universal Lab Analog	13.3mm Length	Titanium	19,-*
C-006-010003		ICX-Lab Analog, short, CAD/CAM	10mm Length	Titanium	19,-*

Note: H = Height from implant shoulder

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)


Overview: ICX-Titanium-Abutments

Article-No.		Article	Dimensions	Material	Price Euro
C-007-010500		ICX-Standard-Titanium-Abutments, Ø 3.7 x 0mm · Not to use with ICX-plus!	0 mm H K 5mm	Titanium	39,-*
C-007-010501		ICX-Titanium-Abutment, short, Standard Ø 5 x 1mm H	1 mm H K 4mm	Titanium	39,-*
C-007-010502		ICX-Titanium-Abutment, short, Standard, Ø 5 x 2mm H	2 mm H K 4mm	Titanium	39,-*
C-007-010503		ICX-Titanium-Abutment, short, Standard, Ø 5 x 3mm H	3 mm H K 4mm	Titanium	39,-*
C-007-010504		ICX-Titanium-Abutment, short, Standard Ø 5 x 4mm H	4 mm H K 4mm	Titanium	39,-*
C-007-010505		ICX-Titanium-Abutment, short, Standard Ø 5 x 5mm H	5 mm H K 4mm	Titanium	39,-*
C-007-010506		ICX-Titanium-Abutment, short, Standard Ø 5 x 6mm H	6 mm H K 4mm	Titanium	39,-*
C-007-016501		ICX-Titanium-Abutment, short, broad Ø 6.5 x 1mm H	1 mm H K 4mm	Titanium	39,-*
C-007-016502		ICX-Titanium-Abutment, short, broad Ø 6.5 x 2mm H	2 mm H K 4mm	Titanium	39,-*
C-007-016503		ICX-Titanium-Abutment, short, broad Ø 6.5 x 3mm H	3 mm H K 4mm	Titanium	39,-*
C-007-020501		ICX-Titanium-Abutment, long, Standard Ø 5 x 1mm H	1 mm H K 7mm	Titanium	39,-*
C-007-020502		ICX-Titanium-Abutment, long, Standard, Ø 5 x 2mm H	2 mm H K 7mm	Titanium	39,-*
C-007-020503		ICX-Titanium-Abutment, long, Standard, Ø 5 x 3mm H	3 mm H K 7mm	Titanium	39,-*
C-007-020504		ICX-Titanium-Abutment, long, Standard Ø 5 x 4mm H	4 mm H K 7mm	Titanium	39,-*
C-007-020505		ICX-Titanium-Abutment, long, Standard Ø 5 x 5mm H	5 mm H K 7mm	Titanium	39,-*
C-007-020506		ICX-Titanium-Abutment, long, Standard Ø 5 x 6mm H	6 mm H K 7mm	Titanium	39,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Note: H = Height from implant shoulder
K = Length of the stump starting from
preparation margin


 All articles on this page: Connection Screw included

 All items on this page: Lab Screw included


Overview: ICX-Titanium-Abutments

Artikel-Nr.		Artikel	Dimensionen	Material	Price Euro
C-007-026501		ICX-Titanium-Abutment, long, broad Ø 6.5 x 1mm H	1 mm H K 7mm	Titanium	39,-*
C-007-026502		ICX-Titanium-Abutment, long, broad Ø 6.5 x 2mm H	2 mm H K 7mm	Titanium	39,-*
C-007-026503		ICX-Titanium-Abutment, long, broad Ø 6.5 x 3mm H	3 mm H K 7mm	Titanium	39,-*
C-007-020500		ICX-Titanium-Abutment, 0mm long, Standard Ø 3.7 x 7mm H Not to use with ICX-plus!	0 mm H K 7mm	Titanium	39,-*
C-007-030500		ICX-Titanium-Abutment, 0mm extralong, Standard, Ø 3.7 x 9mm H Not to use with ICX-plus!	0 mm H K 9mm	Titanium	39,-*
C-010-020015		ICX-Aesthetic Line Titanium Abutment, 15° Aesthetic 1.5mm	1.5mm H	Titanium	79,-*
C-010-020025		ICX-Aesthetic Line Titanium Abutment, 15° Aesthetic 2.5mm	2.5mm H	Titanium	79,-*
C-010-020035		ICX-Aesthetic Line Titanium-Abutment, 15° Aesthetic 3.5mm Titanium	3.5mm H	Titanium	79,-*
C-010-020050		ICX-Aesthetic Line Titanium-Abutment, 15° Aesthetic 5mm Titanium	5.0mm H	Titanium	79,-*
C-010-040015		ICX-Aesthetic Line Titanium-Abutment, 25° Aesthetic 1.5mm Titanium	1.5mm H	Titanium	79,-*
C-010-040025		ICX-Aesthetic Line Titanium-Abutment, 25° Aesthetic 2.5mm Titanium	2.5mm H	Titanium	79,-*
C-010-040035		ICX-Aesthetic Line Titanium-Abutment, 25° Aesthetic 3.5mm Titanium	3.5mm H	Titanium	79,-*
C-010-040050		ICX-Aesthetic Line Titanium-Abutment, 25° Aesthetic 5.0mm Titanium	5.0mm H	Titanium	79,-*
C-010-010010		ICX-Aesthetic Line Titanium Abutment, Aesthetic 1.0mm	1.0mm H K 7mm	Titanium	79,-*
C-010-010030		ICX-Aesthetic Line Titanium Abutment, Aesthetic 3.0mm	3.0mm H K 7mm	Titanium	79,-*

Note: H = Height from implant shoulder
K = Length of the stump starting from preparation margin






 All articles on this page: Connection Screw included

 All items on this page: Lab Screw included

Overview: ICX-Massive- and Universal-Abutments

ICX-Massive-Abutments




Article-No.	Article	Dimensions	Material	Price Euro	
C-021-063401	ICX-Massive-Abutment-Mini, 1mm H, Ø 4.1mm, onepiece	1mm H 5.5mm K	Titanium	39,-*	
C-021-063402	ICX-Massive-Abutment-Mini, 2mm H, Ø 4.1mm, onepiece	2mm H 5.5mm K	Titanium	39,-*	
C-021-063403	ICX-Massive-Abutment-Mini, 3mm H, Ø 4.1mm, onepiece	3mm H 5.5mm K	Titanium	39,-*	
C-021-065001	ICX-Massive-Abutment-Standard, 1mm H, Ø 5mm, onepiece	1mm H 5.5mm K	Titanium	39,-*	
C-021-065002	ICX-Massive-Abutment-Standard, 2mm H, Ø 5mm, onepiece	2mm H 5.5mm K	Titanium	39,-*	
C-021-065003	ICX-Massive-Abutment-Standard, 3mm H, Ø 5mm, onepiece	3mm H 5.5mm K	Titanium	39,-*	
C-021-000002	ICX-Wrench for Massive abutments and straight OCC- & Fangmann-Abutments	-	High-grade steel	49,-*	

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Note: H = Height from implant shoulder
K = Length of the stump starting from
preparation margin

Note: The torque of 25Ncm should not be exceeded. There is a risk that
cone forces may lock it permanently, if screwd with 20 Ncm or more.

ICX-Universal-Abutments






Article-No.		Article	Dimensions	Material	Price Euro
C-009-100001		ICX Universal Abutment, Ø 5.2, H= 9.5 mm incl. connection and lab screw	9.5mm H	Titanium	39,-*
C-009-100002		ICX Universal Abutment, long, Ø 5.2, H= 14.5 mm incl. connection and lab screw	14.5mm H	Titanium	39,-*
C-009-020001		Meso-Telescope-Abutment- Titanium, Ø 7mm, incl. connection and lab screw	Ø 7mm 9mm H	Titanium	39,-*

Note: H = Height from implant shoulder





*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Overview: ICX-Bar-, Gold- and Plastic-Abutments

ICX-Bar-System-Abutments

Article-No.		Article	Dimensions	Material	Price Euro
C-011-010002		ICX-Bar-System-Gold, Ø 4,5mm, UCLA, non hex, gold abutment for bar	9mm HK 6,3mm K 1.36g Gewicht	Gold	109,-*
C-011-030001		ICX-Bar-System-plastic Ø 4.5mm, UCLA, plastic abutment, hex	9mm HK 6.3mm K	Plastic burn out	29,-*
C-011-030002		ICX-Bar-System-plastic Ø 4.5mm, UCLA, plastic abutment, non hex	9mm HK 6.3mm K	Plastic burn out	29,-*
C-011-020002		ICX-Bar-System-Titanium Ø 4.5mm, non hex	9mm HK 6.3mm K	Titanium	39,-*
C-011-020003		ICX-Bar-System-Titanium Ø 4.5mm, long, non hex	13.5mm HK 6.3mm K	Titanium	39,-*


ICX-Gold- and Plastic-Abutments

Article-No.		Article	Dimensions	Material	Price Euro
C-008-010001		ICX-Gold Abutment hex, Ø 5mm, incl. plastic cap, burn out	-	Gold Ceramicor Cendres & Metaux	109,-*
C-008-010002		ICX-Gold Abutment, non hex, Ø 5mm	-	Gold Ceramicor Cendres & Metaux	109,-*
C-008-030001		ICX-Plastic Abutment hex, Ø 5mm	8mm HK 7mm K	Plastic burn out	29,-*
C-008-030002		ICX-Plastic Abutment, non hex, Ø 5mm	8mm HK 7mm K	Plastic burn out	29,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)



 All articles on this page: Connection Screw included

Note: HK = overall height from implant shoulder
K = Length of the stump starting from preparation margin


 All items on this page: Lab Screw included


Overview: ICX-plus-, and CAD/CAM-Adhesive- and Scanabutments

ICX-plus-Abutments



Article-No.		Article	Dimensions	Material	Price Euro
C-024-010501		ICX-plus - Titanium Abutment short, Ø 5, H= 0 mm, K= 5mm incl. connection and lab screw	5mm K	Titanium	39,-*
C-024-020501		ICX-plus - Titanium Abutment long, Ø 5, H= 0 mm, K= 7.5mm incl. connection and lab screw	7.5mm K	Titanium	39,-*

Note: K = Length of the stump starting from preparation margin

 Connection Screw included





 Lab Screw included

ICX-Fangmann-Abutments

Article-No.		Article	Dimensions	Material	Price Euro
C-009-990002		ICX-Fangmann-Abutment, one piece, 5mm	5mm H Ø 6.5	Titanium	39,-*
C-009-990003		ICX-Fangmann-Abutment, one piece, 7mm	7mm H Ø 6.5	Titanium	39,-*

Note: H = Height from implant shoulder


ICX-CAD/CAM-Adhesive- and Scanabutments

Article-No.		Article	Dimensions	Material	Price Euro
C-029-000002		ICX-Adhesive Abutment hex, Ø 4.4mm plus rotation protection	4.8mm H 4.0mm K	Titanium	39,-*
C-029-000004		ICX-Adhesive Abutment hex, GH 0mm, Ø 4.4mm, plus rotation protection · Not to use with ICX-plus!	4.5mm K 4.0mm H	Titanium	39,-*
C-029-001001		ICX-Bar-Adhesive Abutment, Ø 3.8mm	4.3mm K 3.8mm H	Titanium	39,-*
C-030-000001		ICX-Scanabutment	12mm H	Plastic	39,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Note: H = Height from implant shoulder
K = Length of the stump starting from preparation margin

 incl. Connection Screw

 incl. Lab Screw

ICX-Ceramic-Abutments

Article-No.		Article	Dimensions	Material	Price Euro
C-016-020015		ICX-15° Aesthetic Ceramic 1.5mm	1.5mm H	Ceramic	149,-*
C-016-020025		ICX-15° Aesthetic Ceramic 2.5mm	2.5mm H	Ceramic	149,-*
C-016-010501		ICX-Ceramic-Standard-Abut- ment Ø 5mm x 1mm H	1mm H 4mm K	Ceramic	129,-*
C-016-010502		ICX-Ceramic-Standard-Abut- ment Ø 5mm x 2mm H	2mm H 4mm K	Ceramic	129,-*
C-016-010503		ICX-Ceramic-Standard-Abut- ment Ø 5mm x 3mm H	3mm H 4mm K	Ceramic	129,-*

Note: H = Height from implant shoulder
K = Length of the stump starting from
preparation margin






 incl. Connection Screw

 incl. Lab Screw

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)





Overview: ICX-Maximus-Abutments, two- and onepiece

ICX-Maximus-Abutments, twopiece

Article-No.		Article	Dimensions	Material	Price Euro
C-031-010001		1.0 mm ICX-Maximus-Abutment, twopiece incl. Connection screw	1mm Gingiva Height	Titanium	99,-*
C-031-010002		2.0 mm ICX-Maximus-Abutment, twopiece incl. Connection screw	2mm Gingiva Height	Titanium	99,-*
C-031-010003		3.0 mm ICX-Maximus-Abutment, twopiece incl. Connection screw	3mm Gingiva Height	Titanium	99,-*
C-031-010004		4.0 mm ICX-Maximus-Abutment, twopiece incl. Connection screw	4mm Gingiva Height	Titanium	99,-*
C-031-010005		5.0 mm ICX-Maximus-Abutment, twopiece incl. Connection screw	5mm Gingiva Height	Titanium	99,-*

Compatible with LOCATOR-Abutments, Maximus Abutments are an inhouse development of medentis medical

ICX-Maximus-Abutments, onepiece

Article-No.		Article	Dimensions	Material	Price Euro
C-031-020001		ICX-Maximus-Abutment, 1mm, onepiece	1mm Gingiva Height	Titanium	99,-*
C-031-020002		ICX-Maximus-Abutment, 2mm, onepiece	2mm Gingiva Height	Titanium	99,-*
C-031-020003		ICX-Maximus-Abutment, 3mm, onepiece	3mm Gingiva Height	Titanium	99,-*
C-031-020004		ICX-Maximus-Abutment, 4mm, onepiece	4mm Gingiva Height	Titanium	99,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)





Compatible with LOCATOR-Abutments, Maximus Abutments are an inhouse development of medentis medical




Overview: Locator™ for ICX-templant®

Article-No.		Article	Material	Price Euro
ICX-8351		1.0 mm Locator™-Abutment for ICX templant®	Titanium	119,-*
ICX-8352		2.0 mm Locator™-Abutment for ICX templant®	Titanium	119,-*
ICX-8353		3.0 mm Locator™-Abutment for ICX templant®	Titanium	119,-*
ICX-8354		4.0 mm Locator™-Abutment for ICX templant®	Titanium	119,-*
LO-8527		Retention Replacement Male pink (4 pcs.), Retention forces approx. 1.360g, Divergence balance $\pm 10^\circ$	Plastic	19,-*
LO-8547		Retention Replacement Male green (4 pcs.), Retention forces approx. 1.590g, Divergence balance $\pm 20^\circ$	Plastic	19,-*
LO-8524		Retention Replacement Male transparent (4 pcs.), Retention forces approx. 2.265g, Divergence balance $\pm 10^\circ$	Plastic	19,-*
LO-8529		Retention Replacement Male blue (4 pcs.), Retention forces approx. 680g, Divergence balance $\pm 10^\circ$	Plastic	19,-*
LO-8548		Retention Replacement Male red (4 pcs.), Retention forces approx. 680g, Divergence balance $\pm 20^\circ$	Plastic	19,-*
LO-8519-2		Locator™-Male Processing Package (Titanium Abutments with Retention Replacements, 2 pcs.)		43,50*
LO-85-19-10		Locator™-Male Processing Package (Titanium Abutments with Retention Replacements, 10 pcs.)		199,-*
LO-8393		Locator™ Core Tool		66,-*
LO-8505		Locator™-Impression Post (4 pcs.)		59,-*
LO-8530		Locator™-Model Analog (Ø 4 mm) (4 pcs.)		59,-*
LO-8913		Torque Wrench Driver (short) 23 mm		49,-*
LO-8914		Torque Wrench Driver (long) 29 mm		49,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)






Overview: ICX-t-bona Ball-Head-System

Article-No.		Article	Dimensions	Material	Price Euro
C-022-060001		ICX-t-bona Ball-Head-System, 1mm H, ball-Ø 2.25mm incl. Lab Analog	1mm H	Titanium	59,-*
C-022-060002		ICX-t-bona Ball-Head-System, 2mm H, ball-Ø 2.25mm incl. Lab Analog	2mm H	Titanium	59,-*
C-022-060003		ICX-t-bona Ball-Head-System, 3mm H, ball-Ø 2.25mm incl. Lab Analog	3mm H	Titanium	59,-*
C-022-060004		ICX-t-bona Ball-Head-System, 4mm H, ball-Ø 2.25mm incl. Lab Analog	4mm H	Titanium	59,-*

Article-No.		Article	Dimensions	Material	Price Euro
T-13825		T-Ecco, Matrix, complete, including Retention Replacement Male gold and red Plastic Auxiliary Ring	-	Titanium	49,-*
T-13833		T-Ecco, Retention Replacement golden	-	Titanium	35,-*
T-10571		T-Ecco-Aktivator/Deaktivator	-	Titanium	49,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Overview: ICX-mini Ø 2.9mm/Prosthetics

Article-No.		Article	Dimensions	Material	Price Euro
C-026-010501		ICX-mini Massive Abutment (The massive mini abutment is to be cemented on top of the ICX-mini ballhead implant. The cementable mini abutment can be ground afterwards. A conventional impression is made directly from the ground abutment)	Ø 5mm 7mm K	Titanium	39,-*
C-003-000004		ICX-mini Cover Cap, (Note 1!) Only for 2.9mm Implants	-	Titanium	16,50*
C-022-090001		t-bona-Screw, <i>mini</i> , 1mm H, ball-Ø 2.25mm, incl. Model Analog, Only for 2.9mm Implants	1mm H	Titanium	59,-*
C-022-090002		t-bona-Screw, <i>mini</i> , 2mm ball-Ø 2.25mm, incl. Model Analog, Only for 2.9mm Implants	2mm H	Titanium	59,-*
C-022-090003		t-bona-Screw, <i>mini</i> , 3mm ball-Ø 2.25mm, incl. Model Analog, Only for 2.9mm Implants	3mm H	Titanium	59,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Note1: Please tighten all ICX-cover screws and healing caps by hand only. A torque of 5 to 10Ncm should not be exceeded. There is a risk that cone forces may lock it permanently, if screw with 20 Ncm or more.

Note: H = Height from implant shoulder
K = Length of the stump starting from preparation margin

OCC-Product-Line = occlusally screw-retained implant prosthetics

necessary minimum overall height 9.7mm

Article-No.		Article	Dimensions	Material	Price Euro
C-020-150010		OCC Abutment straight 1mm Titanium one-piece, incl. Plastic bringing in assistance	1mm H 4.8mm K	Titanium	39,-*
C-020-150020		OCC Abutment straight 2mm Titanium one-piece, incl. Plastic bringing in assistance	2mm H 4.8mm K	Titanium	39,-*
C-020-150030		OCC Abutment straight 3mm Titanium one-piece, incl. Plastic bringing in assistance	3mm H 4.8mm K	Titanium	39,-*
C-020-152010		OCC Abutment, 20° 1mm Titanium, with hex incl. Lab- and Connection Screw and adjusting assistance	1mm H 4.8mm K	Titanium	99,-*
C-020-153010		OCC Abutment 30° 1mm Titanium, with hex incl. Lab- and Connection Screw and adjusting assistance	1mm H 4.8mm K	Titanium	99,-*
C-020-154010		OCC Abutment 40° 1mm Titanium, with hex incl. Lab- and Connection Screw and adjusting assistance	1mm H 4.8mm K	Titanium	99,-*
C-020-252010		OCC Abutment 20° 1mm Titanium, with hex incl. Lab- and Connection Screw and adjusting assistance	1mm H 4.8mm K	Titanium	99,-*
C-020-253010		OCC Abutment 30° 1mm Titanium, with hex incl. Lab- and Connection Screw and adjusting assistance	1mm H 4.8mm K	Titanium	99,-*
C-020-254010		OCC Abutment 40° 1mm Titanium, with hex incl. Lab- and Connection Screw and adjusting assistance	1mm H 4.8mm K	Titanium	99,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Note: H = Height from implant shoulder

K = Length of the stump starting from preparation margin

 incl. Connection Screw















 incl. Lab Screw

Application and indication:

The OCC-concept is a prosthetic concept for immediate restoration of edentulous patients. For a better support of the provisional restoration it is possible to place the posterior implants with a 20° or 30° angulation to distal. This angulation is compensated by the special OCC-abutments.

OCC-Product-Line = occlusally screw-retained implant prosthetics




















necessary minimum overall height 9.7mm

Article-No.		Article	Dimensions	Material	Price Euro
C-020-351125		OCC Prosthetic Abutment, 12.5mm, Titanium, incl. Connection Screw	12.5mm K	Titanium	39,-*
C-020-352125		OCC Prosthetic Abutment, 12.5mm, plastic burn out, incl. Connection Screw	12.5mm K	plastic burn out	19,-*
C-020-353125		OCC Prosthetic Abutment, gold incl. burn outable plastic cap, incl. Connection Screw	12.5mm K	Gold/HSL-Alloy Ceramicor Cendres & Metaux	109,-*
C-003-000003		OCC Protection cap, Ø 5mm, incl. Connection Screw	6.6mm K	Titanium	19,-*
C-020-381125		OCC Bar-Abutment, Ø 5mm Titanium, incl. Connection Screw	8.8mm K	Titanium	39,-*
C-020-382125		OCC Bar-Abutment, Ø 5mm plastic, incl. Connection Screw	8.8mm K	Plastic burnout	19,-*
C-020-392125		OCC Interim/Bite Registration Peek, incl. Connection Screw	12.5mm K	Peek	19,-*
C-005-050001		OCC Impression Post closed, incl. Casting Cap incl. Connection Screw	Article replaces the past article with the same article number!	Titanium	29,-*
C-005-050002		OCC Impression Post open, incl. Connection Screw	7.7mm K	Titanium	29,-*
C-006-010004		OCC Model Analog straight, Ø 5mm	17mm Length	Titanium	19,-*
C-006-010006		OCC Model Analog short, CAD/CAM, Ø 5mm	10.8mm length	Titanium	19,-*
C-015-100016		OCC Angle Template for 30° and 40°	-	Titanium	79,-*
C-029-002001		OCC Adhesive Base, incl. Connection Screw	9mm K	Titanium	39,-*
C-030-000002		OCC Scanbase, incl. Connection Screw	12.5mm K	plastic	39,-*

Note: H = Height from implant shoulder
K = Length of the stump starting from preparation margin

 incl. Connection Screw

Overview: ICX-Parallel- and Stop-Drills

Article-No.	Drills	Article	Price Euro
FIL-186RF		ICX-Cortical-Drill	59,-*
C-014-006290		ICX-Parallel-Drill WHITE = Ø 2.9mm	59,-*
C-014-007290		ICX-Parallel-Drill WHITE = Ø 2.9mm	59,-*
C-014-003345		ICX-Parallel-Drill ● = Ø 3.45mm	59,-*
C-014-006345		ICX-Parallel-Drill ● = Ø 3.45mm	59,-*
C-014-007345		ICX-Parallel-Drill ● = Ø 3.45mm	59,-*
C-014-003375		ICX-Parallel-Drill ● = Ø 3.75mm	59,-*
C-014-006375		ICX-Parallel-Drill ● = Ø 3.75mm	59,-*
C-014-007375		ICX-Parallel-Drill ● = Ø 3.75mm	59,-*
C-014-003410		ICX-Parallel-Drill ● = Ø 4.1mm	59,-*
C-014-006410		ICX-Parallel-Drill ● = Ø 4.1mm	59,-*
C-014-007410		ICX-Parallel-Drill ● = Ø 4.1mm	59,-*
C-014-003480		ICX-Parallel-Drill ● = Ø 4.8mm	59,-*
C-014-006480		ICX-Parallel-Drill ● = Ø 4.8mm	59,-*
C-014-007480		ICX-Parallel-Drill ● = Ø 4.8mm	59,-*
C-014-005345		ICX-Stop-Drill ● = Ø 3.45mm	59,-*
C-014-005375		ICX-Stop-Drill ● = Ø 3.75mm	59,-*
C-014-005410		ICX-Stop-Drill ● = Ø 4.1mm	59,-*
C-014-005480		ICX-Stop-Drill ● = Ø 4.8mm	59,-*
C-014-000003		ICX-Crestotom-Drill	59,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)



59,€*
per ICX-implant
all lengths
all diameters
*plus VAT

**SILVERSPONSOR
of DGI**

**BRONZESPONSOR
of EAO**

Overview: ICX-Drillstops, Titanium

Article-No.	Article	Material	Price Euro	
C-027-129080	ICX-Drillstop Ø 2.9mm, Titanium, Drill depth 8mm Colour code: ○ WHITE	Titanium	19,90*	
C-027-129100	ICX-Drillstop Ø 2.9mm, Titanium, Drill depth 10mm Colour code: ○ WHITE	Titanium	19,90*	
C-027-129125	ICX-Drillstop Ø 2.9mm, Titanium, Drill depth 12.5mm Colour code: ○ WHITE	Titanium	19,90*	
C-027-135080	ICX-Drillstop Ø 3.45mm, Titanium, Drill depth 8mm Colour code: ● YELLOW	Titanium	19,90*	
C-027-135100	ICX-Drillstop Ø 3.45mm, Titanium, Drill depth 10mm Colour code: ● YELLOW	Titanium	19,90*	
C-027-135125	ICX-Drillstop Ø 3.45mm, Titanium, Drill depth 12.5mm Colour code: ● YELLOW	Titanium	19,90*	
C-027-138080	ICX-Drillstop Ø 3,75mm, Titanium, Drill depth 8mm Colour code: ● RED	Titanium	19,90*	
C-027-138100	ICX-Drillstop Ø 3,75mm, Titanium, Drill depth 10mm Colour code: ● RED	Titanium	19,90*	
C-027-138125	ICX-Drillstop Ø 3,75mm, Titanium, Drill depth 12.5mm Colour code: ● RED	Titanium	19,90*	
C-027-141080	ICX-Drillstop Ø 4,1mm, Titanium, Drill depth 8mm Colour code: ● GREEN	Titanium	19,90*	
C-027-141100	ICX-Drillstop Ø 4,1mm, Titanium, Drill depth 10mm Colour code: ● GREEN	Titanium	19,90*	
C-027-141125	ICX-Drillstop Ø 4,1mm, Titanium, Drill depth 12.5mm Colour code: ● GREEN	Titanium	19,90*	
C-027-148080	ICX-Drillstop Ø 4,8mm, Titanium, Drill depth 8mm Colour code: ● BLUE	Titanium	19,90*	
C-027-148100	ICX-Drillstop Ø 4,8mm, Titanium, Drill depth 10mm Colour code: ● BLUE	Titanium	19,90*	
C-027-148125	ICX-Drillstop Ø 4,8mm, Titanium, Drill depth 12.5mm Colour code: ● BLUE	Titanium	19,90*	















*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Article-No.	Article	Material	Price Euro
C-027-000010	ICX-Drillstop-Box (empty)	Plastic, sterilizable	99,-*
C-027-000100	ICX-Drillstop-Box (incl. 15 Drillstops-Titanium)	Plastic, sterilizable	349,-*










*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)






Overview: ICX-templant®-Accessoires

Article-No.		Article	Material	Price Euro
C-015-100001		ICX Insertion Tool, short, Length 23mm Insertion Tool short with ISO-shaft For all Ø 3.45/3.75/4.1/4.8er implants	High-Grade-Steel	29,-*
C-015-100002		ICX Insertion Tool, medium, Length 29mm Insertion tool medium-sized with with ISO-shaft, For all Ø 3.45/3.75/4.1/4.8er implants	High-Grade-Steel	29,-*
C-015-100004		ICX Insertion Tool long, Length 40mm Insertion tool long with with ISO-shaft For all Ø 3.45/3.75/4.1/4.8er implants	High-Grade-Steel	29,-*
C-015-100003		ICX Depth Measuring Tool For inspection of the drilling depth 6.5/8/10/12.5/15mm	Titanium	79,-*
C-015-100005		ICX Extraction Tool	High-Grade-Steel	49,-*
C-015-100006		ICX Insertion Tool, ISO-Shaft with 1mm T-Bona screw exclusively (standard)	High-Grade-Steel	49,-*
C-015-100007		ICX-t-bona Insertion Tool (ISO-shaft) Installation of all t-Bona screws	High-Grade-Steel	49,-*
960001		Ratchet	High-Grade-Steel	159,-*
960004		Universal Drill Extension for all ICX-Drills	High-Grade-Steel	49,-*
960007		Universal Manual Adapter	High-Grade-Steel	49,-*
950099		Hex Tool 1.4 medium, Length 24mm for all Abutment Screws	High-Grade-Steel	29,-*
950098		Hex Tool 1.4 long, Length 28mm for all Abutment Screws	High-Grade-Steel	29,-*
950097		Hex Tool 1.4 extra long, Length 45mm, for all Abutment Screws	High-Grade-Steel	29,-*
C-021-000002		ICX-Wrench for Massive Abutments and straight OCC- & Fangmann-Abutments	High-Grade-Steel	49,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

Overview: ICX-templant®-Accessoires

Article-No.		Article	Material	Price Euro
960002		Hex Tool 0.9 for T-ecco female inner part and ICX-mini cover screw	High-Grade- Steel	29,-*
C-015-100021		ICX Insertion Tool, medium, Insertion Tool medium-sized with with ISO- shaft, For all Ø 3.45/3.75/4.1/4.8er implants	High-Grade- Steel	49,-*
C-015-100025		Hex Tool 1.4 medium, for all Abutment Screws	High-Grade- Steel	49,-*
C-015-100020		ICX-Surgical Driver, for instruments with dental shaft	High-Grade- Steel	99,-*
TW-1000003-10		ICX Abutment Mounting, threadable into handhold TW-100100	High-Grade- Steel	49,-*
TW-1000100		Handhold, length 130mm, Ø 15mm Laboratory Tool for Abutment grinding (without abutment mounting)	High-Grade- Steel	49,-*
C-015-100009		Countertool for <i>ICX-mini</i> (for ballhead removal)	High-Grade- Steel	29,-*

Article-No.		Article	Price Euro
C-014-230315		ICX Parallelpin Ø 2.3/Ø 3.15 for drilling cavity control	29,-*
C-014-350410		ICX Parallelpin Ø 3.5/Ø 4.1 for drilling cavity control	29,-*
C-014-006002		ICX Boneprofiler for excavation of covered implants	59,-*
960004		Universal Drill Extension for all ICX-Drills	49,-*
C-015-100003		ICX Depth Measuring Tool For inspection of the drilling depth 6.5/8/10/12.5/15mm	79,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)

ICX-Surgery boxes

Article-No.	Article	Price Euro
-	<p>ICX-templant®-STARTERSET <i>MINI</i> (54 ICX-templant® x € 59,- = € 3.186 + € 704,- Special price for Surgical-Box) You will receive 54 templants® (please mark the sizes needed), you will also receive a complete surgical set consisting of full set of drills, torque wrench, surgical box and all hand instruments (3 Sets) on permanent loan. We ask for recall for the selection of the implants. This supply is valid only for new customers, up to revocation.</p>	3.890,- *
-	<p>ICX-templant®-STARTERSET <i>PROFESSIONAL</i> You will receive 84 templants® (please mark the sizes needed), you will also receive a complete surgical set consisting of full set of drills, torque wrench, surgical box and all hand instruments (3 Sets) on permanent loan. We ask for recall for the selection of the implants. This supply is valid only for new customers, up to revocation.</p>	4.956,-*

*All prices in Euro € plus VAT, until further notice (Date June, 01. 2012)



ICX *templant*[®]
www.medentis.de



59,€*
per ICX-implant
all lengths
all diameters
*plus VAT

www.medentis.de

ICX *templant*[®]
www.medentis.de

ICX-Systemoverview 6.0

... reloaded ...



59,€*
per ICX-implant
all lengths
all diameters
*plus VAT

www.medentis.de